

CRITICAL CARE/PULMONARY THERAPIES

Critical Care Patient

Pulmonary issues

Bariatric Patient >300 lbs

ARDS

Yes

No

Unstable cervical/
thoracic/lumbar fx or
skeletal traction

Yes



RotoRest
Max weight 350 lbs

Indications:

- Treatment and prevention of pulmonary complications due to immobility
- Treatment of ARDS



Triadyne II
Max weight 300 lbs

Indications:

- Treatment and prevention of pulmonary complications associated with immobility
- *Order Triadyne Proventa for patients up to 300 lbs.**



BariAir
Max weight 850 lbs

Indications:

- Low air loss/turn assist
- Patient whose weight poses a risk to patient and staff
- Preventing and treating pressure ulcers
- Large patients requiring percussion therapy
- ICU ONLY**



Roto-Rest Delta
Max weight 300

Indications:

- Treatment and prevention of pulmonary complications due to immobility
- Thoracic and lumbar fracture
- Cervical traction (not as a primary means)
- Skeletal traction